


LCMS Identification using LipidSearch 5.0 Documentation

LipidSearch is designed to process LCMS data and perform lipid identification using data-dependent MS2 data.

1.1 Product Search(LCMS) Overview

LipidSearch 5.0 performs lipid identification using the following procedure.

1. Processes the Extracted Ion Chromatogram [EIC] and detects peaks [Separated Peak Area: SPA] for each lipid ion.
2. Collects the MS2 spectra contained in the EIC region.
3. Calculates the Product Ion Chromatogram [PIC] of the product ion peak in each MS2 spectrum and calculates the virtual spectrum [Deconvoluted MS2 Spectrum: DMS] distributed from the shape belonging to the SPA.
4. Lists lipid ions with matching DMS and peak patterns as candidates.



1.2 Identification Metrics

Two types of indicators are used to ensure the certainty of identification.

1.2.1 Grade

Grade is calculated as follows using the product ion Category specified in the Database.

- [Category]

Category	Description
C0	Head group specific ions

Category	Description
C01	Head group-derived or substituent-derived ions specific to lipid class
C1	Substituent-specific ions
C2	Ions that do not identify structure

- [Condition]

Condition	Description
H1	C0 and C01 (or more than one C0) (or C01 for all substituents) are assigned
H2	C0 or C1 are assigned
S1	C01 or C1 are assigned for all substituents
S2	C01 or C1 are assigned for some substituents

- [Grade Assignment]

Grade	Description
A	(H1 and S1)
B	(H2 and S1) or (H1 and S2)
C	H2 or S1
D	Other

1.2.2 IDScore

The identification score indicates to what degree the product ion peaks are assigned to a given identification. The more a peak of high intensity in [DMS](#) is specifically associated with lipid candidates, the higher the score. If all the peaks in DMS are assigned to only one lipid candidate, the score is 1.0, and if there are many peaks that are assigned across multiple candidates, the score is lower. In some cases, if there are many peaks that are only assigned for one lipid among multiple candidates, it exceeds 1.0.

$$IDScore = \frac{1}{N_s} \sum_k \left(\left(\frac{I_{m(i,k)}}{\sum_j I_{(i,j)}} \right) \cdot \log \left(\frac{N_l}{Lf(P_k)} + 1 \right) \right)$$

$I_{m(i,k)}$: Intensity of peak k consistent with product ions of lipid candidates in spectrum i

$I_{(i,j)}$: Intensity of peak j in spectrum i

N_s : Number of spectra used to identify these lipid candidates (1+MS3 spectra)

$Lf(p_k)$: Number of lipid candidates with product ions assignment to peak k

N_l : Number of lipid candidates for this DMS

Here are the parameters used for identification

Parameter

Identification | Database Setting | Manual Filter | Quantitation

☐ Save options DE12-3-POS

▼ Identification

SearchType Product Search
ExpType LC-MS
Precursor tol 5.0 ppm
Product tol 5.0 ppm
Product Threshold 1.0 relative

▼ ID Range

☒ rt 1.0 ~ 26.0
☒ precursor m/z 300.0 ~ 1800.0

▼ Others

Peak Detection

Peak Detection Parameter Settings

BUNCHING_FACTOR.MS1	5
NOISE_STD.PER.MS1	0.001
NOISE_STD.TIME.MS1	0.01
BUNCHING_FACTOR.MSN	3
NOISE_STD.PER.MSN	0.0
NOISE_STD.TIME.MSN	0.01
ABSL_INTENSITY_THRESHOLD.MS1	30000.0
SECDENUM_RATIO_THRESHOLD.TO.MAX	0.3
SAME_PEAK_RT_THRESHOLD	0.025
QUANT_CLUSTER_COS_THRESHOLD	0.9
QUANT_INTENSITY_RATIO_THRESHOLD	3.0
CONNECTIVITY_RATIO_THRESHOLD	0.05
CHARGE_DETECT_DELTA.MZ	0.02
MIN_PEAK_CHARGE	1
MAX_PEAK_CHARGE	2
CHARGE_DETECT_COS_THRESHOLD	0.9
CHARGE_DETECT_RECALC_THRESHOLD	0.7
CHROM_RET_INTENSITY_THRESHOLD	0.005
DO_BASLINE_SUBTRACT	0
BASLINE_TO_MAX_RATIO_THRESHOLD	1.0
CHROM_FITTING_GAUSS	0
CHARGE_DETECT_MIN_PEAK_NUM	2